

Table of Contents

Page

In vivo stability of WS1321 ASD strains carrying HIV peptides	1-2
Cloning and testing HIV genes under the control of pagC	3-4
Macrophage Assay of pagC Construct	5
Growth of WS1321 ASD + pAB104 in Minimal Media Supplemented w/ 0.4% glucose	6
Mouse Experiment pagC constructs	7
Construction of WS1321 ASD strains carrying pAB105 and pAB105::vif for mice feeding	8
Construction of pAB105::ENV110	9
Summary of Human Monocyte & Staphy Interaction Exp Trial (P78)	10
Cloning of ENV110-A1 (cont.)	11-12
Shigella CMI Response	13
Shigella as a DNA Carrier (Antibiotic Concept)	14
Shigella as a DNA Carrier (Burno Jirlands Concept)	15
Shigella as a DNA Carrier (Vectors based on RNA)	16
" - Monkey Experiments	17
Use of amphotericin to Kill Attenuated Shigella Strains & Allow Delivery of Plasmid DNA to Eukaryotic Cells	18
Antibiotic Use to Kill Attenuated Shigella Strain SG602 in Order to Facilitate Delivery of Plasmid DNA	19-20
	21

Title: In vivo stability of WS1321 AASD strains carrying HIV peptide

Date:

Purpose: Determine the number of bacteria recovered from the spleens of orally infected BALB/c mice. Compare the levels of each ^{mutant} strain to that of the parent strain.

Strains: Constructed by Art Brantstrom are:

WS1321AASD, PAB102

WS1321AASD, PAB103

WS1321AASD, PAB02::gag

WS1321AASD, PAB03::gag

WS1321AASD, PAB02::vif

WS1321AASD, PAB03::vif

Also

WS1321 - plasmid - Sal. typhimurium strain

See to explain purpose of construction of these 2 strains.

Note: A previous experiment conducted indicated by day 6 after ^{oral} infection fewer numbers of strain WS1321, PAB103 were recovered from the spleens. ← Lab Notebook #3 (House work)

Experiment: (Full procedure documented in Lab Notebook #3 (House work))

- strains were grown to an O.D. of $\approx 0.8 - 1.0$
- ~~cells~~ to concentrate into DPBS (Some cultures resuspended in 0.35ml. others into 0.5ml)
- cultures were diluted & plated to determine ^{amt of} bacteria fed.
- 40ul fed to mice using 200ul pipet.
- mice are to be sacrificed on Days 4, 7 & 11 post-feeding
- 3 mice each strain each timept.

Amt. Given to Mice: Between $3 - 4 \times 10^9$ in 40ul.

Exact #'s are in Book #3

* Spleens homogenized in 2ml HBSS 0.1ml original plated in duplicate

Continued on Page 2

Amata R. Lyons

Signed

Date

Rear REDACTED

Signed

Date

PROJECT Salmonella Va ie Canis

Continued From Page 1

Results:		ave/splen = ave of 3 x 1/101 x 2nd dil.				
Strain	Day	total ave/splen	Day	total ave/spl.	Day	total ave/spl.
WS1321 #1	4	90, 101	7	128, 122	11	98, 82
#2	131, 46	> 1697	98, 77	> 2097	111, 103	> 1690
#3	70, 71		107, 97		65, 51	
WS1321 #1	44, 37	215	8, 9		63, 100	
PAB102 #2	14, 86	> 1160	92, 135	> 1317	19, 17	> 1123
#3	18, 66		73, 78		74, 67	
WS1321 PAB103 #1	5, 4		35, 46		97, 57	
#2	2, 9	> 70	59, 45	> 767	2, 3	> 1283
#3	0, 0		26, 21		104, 122	
WS1321 PAB102 #1	0, 0		0, 0		2, 4	
#2	0, 0	> 0	18, 19	> 507	0, 0	> 46
#3	0, 0		54, 61		5, 1	
WS1321 #1	0, 0		0, 0		0, 0	
#2	0, 0	> 0	0, 0	> 0	3, 0	> 10
#3	0, 0		0, 0		0, 0	
WS1321 PAB102 #1	0, 0		0, 0		1, 2	
#2	0, 0	> 0	1, 1	> 27	6, 9	> 63
#3	0, 0		0, 6		4, 0	
WS1321 PAB103 #1	0, 0		0, 0		1, 0	
#2	0, 0	> 0	0, 0	> 0	4, 4	> 53
#3	0, 0		0, 0		3, 4	
> 20 not detectable						

Continued on Page

Continued on Page X

Read and Understood By

Signed

Date

Signed

Date

Cloning and Testing HIV genes under the control of the
pagC promoter

Art Bianchini constructed plasmid pAB104, which
is plasmid pAB103 with the pagC promoter replacing
the lac promoter.

* All data in Notebook #1 Digestions & Ligations
Using directed cloning techniques a 2.7 Kb fragment
encoding VIF from HIV-1 was cloned into pAB104

Buildup:

pAB104 2ul
B buffer 2ul
BamHI 1ul
dH₂O 1ul
20ul

pAB104::VIF 10ul
B buffer 2ul
BamHI 1ul
dH₂O 7ul
20ul

↓ 1 hr. 37°C

- add to each & continue at 37°C
for several hours.

above 20ul
H buffer 3.2ul
SalI 2ul
dH₂O 14.8ul
40ul

- + Each rxn was run on a 1.2% low melt agarose
gel. 1/2 rxn mix on each side of the gel. 1/5th
gel was stained w EtBr after electrophoresis. The
appropriate bands were sliced from the unstained side
- Each band was Prep-A-Mix purified & run
through a Millipore column for purification &
resuspended in 25ul TE for ligation

Continued on Page

4

Renato R. Lugo
Signed

Read and Understood By

REDACTED

Date

Signed

Date

Ligation:

Control * PAB104 cut, purified & ligated

10x buffer 5ul

10mM ATP 4ul

PAB104 cut, purified 12.5ul

T4 ligase 2ul

dH₂O 26.5ul

50ul

PAB104 + vif

10x buffer 5ul

10mM ATP 4ul

PAB104 12.5ul

vif 25ul

T4 Ligase 2ul

dH₂O 15ul

50ul

↓ 15°C, overnight in PCR machine

1/2 mix 1/2 mix gel

- Ethanol prec. & electroporate into X6097 ASD

Time (msec)

colonies

Ligation Mix

PAB104 13.7


0

PAB104::vif 14.8

many

5 positives were selected for digestion to look for inserts:

- Wizard Miniprep - Purify DNA
- Cut w/ BamH - SmaI
- Run a gel.



all five appear to have inserts.

all five were prepared for Western Blot examination of expression levels. They were electroporated into JE501110:ASD for

placement into 1321:ASD

Continued on Page

5

Read and Understood By

Signed

Date

Signed

Date


PROJECT

Salmonella Vaccine Carriers

Continued From Page

Title: Macrophage Assay of pagC construct
 Purpose: Determine whether pagC promotes increases expression in murine Macrophages

Notebook page #1 full experimental protocol.
 Assay based on protocol of O'Brien et al.
 CFU/Macrophage & Expression via Western blot were assayed.
 Results:



Western: in pagC notebook - Shows expression from pagC (pAB104::vif only)
 Slides of Macrophages infected with strains - notebook - pagC #1
 Experiment Repeated.

Continued on Page

Read and Understood

REDACTED

Signed

Date

Signed

Date

Title: Growth of WS1321 ΔASD + pAB104 in Minimal Media
Supplemented w 0.1% Glucose

Notebook page #1

Repeated

Protocol:

- Inoculate colony of each strain into 5ml of LB or MEM.
- Grow 3rd & 4th hrs.
- Concentrate & resuspend in dissociation buffer, vortex freeze & thaw
- Determine protein concentration.
- Run Western: compare LB & MEM expression of plac & pagC promoters.

Results: (Notebook page #1)

~~3/16~~ EXP.

Bacterial Counts also determined

Lane # 10. Marker

LB { 9. WS1321
8 pAB103::vif
7 pAB104
6 pAB104::vif

5 WS1321
4 pAB103::vif
3 pAB104
2 pAB104::vif } MEM

1 2 3 4 5 6 7 8 9 10

1 2 3 4 5 6 7 8 9 10

Read and Understood By

REDACTED

Signed

Date

Signature

D:

PROJECT _____

Continued From Page _____

Title Mouse Experiment page Constructs

Strains Tested: WS1321

WS1321 Δ ASD+ pAB103::vifWS1321 Δ ASD+ pAB104WS1321 Δ ASD+ pAB104::vif

Purpose: Determine the number of bacterial colonies surviving mouse passage & expression of vif gene after passage.
 Mice sacrificed on days 3, 7 & 10

Amount of Bacteria Fed to Mice

40ul

WS1321 7.4×10^9 pAB103::vif 4.2×10^9 pAB104 6.9×10^9 pAB104::vif 6.4×10^9

Results: Experimental problem: Apparently mice were mixed by animal handlers. " " ← Based on Western blot analysis.

Donata L. Azemir

Continued on Page _____

REDACTED

Signed _____

Date _____

Signed _____

Date _____


Title: Construction of WS1321ΔASD strains carrying pAB105 & pAB105::vif for mice feeding

Purpose: pAB105 & pAB105::vif plasmids were constructed by OutBranston. In this plasmid the asd gene is from *E. coli*. We believe strains carrying plasmids to this asd gene are healthier.

5/12 - Mice fed the following strains & amount
46ul

WS1321ΔASD + pAB105	3.6×10^9
WS1321ΔASD + pAB105::vif	2.72×10^9
WS1321ΔASD + pAB102::vif	6.4×10^9

Results: CFU/spleen & expression of vif from recovered colonies



Results of Western
Only pAB102::vif
showed significant
expression.

Continued on Page

Read and Understood

REDACTED

Signed

Date

Signed

Date

Title: Construction of pAB105::ENVITB

Purpose: Of the HIV genes we work with only ENV has mapped CTC epitopes for BALB/c m.c.

After several attempts: → it was determined that the entire ENV gene isn't stable in our pAB105 plasmid. Various deletions resulted after each cloning attempt. Ended July

Will try smaller fragment which has been constructed by Art Blanton

End *Donata H. Jensen*

Continued on Page _____

REDACTED

Understood By _____

Signed _____

Date _____

Date _____

Title: Summary of Human Monocyte & S. typhi
Interaction Experiments. July

Note: This work has been in progress for some time.
This entry is a summary of the work completed
& a list of work needed to be done for a publication.

Continued on Page _____

Read and Understood By _____

Signed _____

Date _____

Signed _____

Date _____

PROJECT

Continued From Page

Trial CTL Assay

Purpose: To standardize non-radioactive CTL assays from Biomega. Testing various Fetal Calf Serum lots.

Testing:

#1 FBS certified (low LDH ≈ 470 units) Gibco
Lot # 40K0242

#2 Hyclone

5ml of #1 or #2
45ml of RPMI w/out Phenol red.
50 ml

- add 50ul to each well
- add 50ul of Control + LDH to 6 wells
- add 10ul of Cyto solution
- add 50ul of Substrate Mix (12ml of Assay Buffer to bottle of Substrate Mix.) Cover plate w foil & incubate at RT for 30'
- Add 50ul of Stop Solution
- Remove large bubbles. Read Ab at 490.

Note - better control RPMI 1640 w/out serum.

Results: #1 best #2 - cannot be used - completely red.
Aug 16 1994

Control Testing Maximum & Spontaneous Release of P815

#1 Group - RPMI

#2 Group - Media #2 - CTL Assay Media

Note: Cells are healthy. This is a must.

5×10^3 & 1×10^4 cells were tested for spontaneous & maximum release.

RPMI w/out phenol red & Media #2 w/out phenol red. were used as control backgrounds.

Continued on Page

12

Read and Understood By

REDACTED

Signed

Date

Date

Results:

	RPMI	Max	Min		Media	+2	
5×10^3	0.479	0.138	$\times 3.5$		0.537	0.153	$\times 3.5$
1×10^4	0.748	0.158	$\times 4.7$		1.083	0.150	$\times 7.22$

$$RPMI = 0.074$$

$$Media = 0.162$$

$\rightarrow 3\% \text{ FBS}$

$$\text{Gibco } 10\% = 0.326$$

Assay looks like it should work.

Trial run for CTL Proliferation & Cytokine Assay in *Shigella* 2a, ENVIS-A experiment.

Animals inoculated by Corey Mallet.

Above data contained in *Shigella* Exp/Hela/Hep-2/Notes.

Based on info by Art Branstetter a pUC18 plasmid encoding ENVIS-A fragment is stably maintained in *Shigella flexneri* 2a without selection.

Art Branstetter has checked all strains used for this study for Congo Red binding, Toxicity & plasmid content.

Continued on Page

Read at REDACTED

Signed

Date

Signed

Date

Cloning of ENVTRA-A₁ into pAB105

Purpose: This ENVTRA-A₁ fragment contains the relevant V3 loop, CTL epitope & major Ab epitope. Furthermore, it appears to be more stable in Euc18.

Outline of cloning procedure in rough notes: Page C Notebook

Results: Appears to be cloned. Several weeks of problems to ligase!

Several Ab used for detection of expression

1. human sera (Ab's)
2. human sera (TD)
3. MAb R 9-2
4. V3 (HIV-V3-13) - only one that is good!

Gels Run

Comparing pAB102::ENVTRA-A₁ (plac) to pAB105::ENVTRA-A₁ (page) +

Repeated Results: - plac lysates from Art for 11-28-94 not enough protein loaded
- page expression is clearly higher than plac.

Need to repeat w/ fresh Ab, but all looks good.
Set-up for mouse testing.

Continued on Page

Read and Understood By

REDACTED

Signed

Date

Date

PROJECT _____

Shigella CMI Response50 mice total: 4th floor

Purpose: Test for CTL activity to EN1 epitope.

Mice sacrificed Sept. 6: Note: animals look ruffled.
- 2 *S. flexneri* 2a died.

Cells prepared for CTL, proliferation & cytotoxic assays.

Concept dropped

Continued on Page

Read and Understood By

REDACTED

Signed _____

Date _____

Signed _____

Previous background experiments appear in
Notebooks 1 & 2.

As of _____, to extend this concept we
will try antibiotic treatment of BHK cells which
have been infected w/ SC602(pCMVp). The
idea is that antibiotic treatment will kill
the intracellular bacteria. The dead bacteria
will release the pCMVp plasmid into the cytoplasm
for transcription & translation by the BHK cell.

This entry is based on a conversation with
Jung Hwang on _____.

- This entry was

Wendy R. Sugrue

Continued on Page

Read and Understood By

Lillian Van O'Vey

Signed

Date

Wendy R. Sugrue

Signed

as seen

Date